






















Codice	00500.00
Descrizione	Kit elettrodi assortiti (10 pz) 5cm
Lunghezza elettrodi	5 cm
Attacco Manipolo	Ø 2,4 mm
1 x L1	 Elettrodo dritto a filo fine
1 x L2	 Elettrodo angolato a filo fine
1 x L3	 Elettrodo ad ansa Ø 4 mm
1 x L4	 Elettrodo ad ansa Ø 8 mm
1 x L5	 Elettrodo angolato ad uncino
1 x L6	 Elettrodo angolato a filo spesso
1 x L7	 Elettrodo a goccia
1 x L8	 Elettrodo a cappio
1 x L9	 Elettrodo dritto a sfera
1 x L10	 Elettrodo angolato a sfera
Sterilizzazione	Autoclave 134 °C (10 minuti) / 121 °C (20 minuti)
Cicli di sterilizzazione	30



Item	00500.00	
Description	Kit of assorted electrode length 5cm (10 pcs)	
Electrode length	5 cm (1,96 in)	
Attack handpiece	Ø 2,4 mm (0,094 in)	
1 x L1		Straight thin wire electrode
1 x L2		Bent thin wire electrode
1 x L3		Loop electrode Ø 4 mm
1 x L4		Loop electrode Ø 8 mm
1 x L5		Hook electrode
1 x L6		Bent thick wire electrode
1 x L7		Drop electrode
1 x L8		Noose electrode
1 x L9		Straight ball electrode
1 x L10		Bent ball electrode
Sterilization	Autoclave 134 °C [275 °F] (10 minutes) / 121 °C [250 °F] (20 minutes)	
Sterilization cycle	30	

Data, descriptions, dimensions, weights and illustrations are indicative and could be modify without warning.